

|  |                              |
|--|------------------------------|
| <b>Technical Specifications Curatron</b>         | <b>2000PC</b>                |
| Pad 50 x 70 cm = 20" x 28"                       | included                     |
| High energy coil. Diameter 20 x 20 cm = 8"       | option                       |
| Very high energy coil. Diameter 20 x 20 cm = 8"  | option                       |
| Mattress 65 x 180 cm = 26" x 71"                 | option                       |
| Max energy Pad                                   | 40 milli-Tesla (400 Gauss)   |
| Max energy high energy coil                      | 70 milli-Tesla (700 Gauss)   |
| Max energy very high energy coil                 | 100 milli-Tesla (1000 Gauss) |
| Max energy mattress                              | 15 milli-Tesla (150 Gauss)   |
| Master programs                                  | 10 + 70                      |
| Program time for 10 programs                     | 20 – 30 min.                 |
| Frequency range                                  | 1 – 50 Hz                    |
| User programmable master programs                | Yes                          |
| Patient database and Therapy protocols databases | Yes                          |
| Minimum stage time for each parameter setting    | 1 min.                       |
| Maximum stage time for each parameter setting    | 20 min.                      |
| Minimum total therapy time per session           | 1 min.                       |
| Maximum total therapy time per session           | Unlimited                    |