

<b>Technical Specifications Curatron</b>	<b>2000XP</b>
Pad 50 x 70 cm = 20" x 28"	included
High energy coil. Diameter 20 x 20 cm = 8"	option
Very high energy coil. Diameter 20 x 20 cm = 8"	option
Mattress 65 x 180 cm = 26" x 71"	option
Max energy Pad	30 milli-Tesla (300 Gauss)
Max energy high energy coil	45 milli-Tesla (450 Gauss)
Max energy very high energy coil	70 milli-Tesla (700 Gauss)
Max energy mattress	10 milli-Tesla (100 Gauss)
Master programs	10
Program time for 10 programs	20 – 30 min.
Frequency range	1 – 50 Hz

Energies over complete coil area **simultaneously**