



Curatronic Ltd.

Health Care Practitioners

Your Partner in Healing



Curatronic 2000PC

With Pre-programmed Treatment Software



**ENERGY
MEDICINE**

Curatron 2000PC

The Curatron 2000PC clinical system combines the professional XP model with a powerful software package for a PC or MAC. This advanced system offers a unique solution for a Comprehensive Clinical Therapy Center and allows for the highest possible energy output of up-to 100,000 microTesla with a frequency range between 1 and 50Hz.

No special training is required for operating this user-friendly system. The strengths of the applied pulsed electromagnetic fields, the unique frequency modulation and the fast speed of induction are uncontested by any other clinical PEMF system. This makes the efficacy of the Curatron 2000PC pulsed magnetic field therapy system far superior!



Targeted Treatment

The Curatron 2000PC model allows the PEMF therapist to easily target specific conditions with the use of the pre-programmed disorder settings (see back page). The energy coils or probes can be used to deliver a focused energy field for localized treatment.

Bone and Cartilage Regeneration

Improved calcium transport increases the absorption of calcium in bones and improves the quality of cartilage in joints, decreasing pain dramatically. Substantial increase in bone density is obtained in cases of osteoporosis, which in turn decreases fracture risk. Bone growth

is stimulated causing fractured bones to heal properly. On the X-ray on the left, the distance between the bones is almost non-existent. After PEMF treatment, the X-ray on the right shows the distance between the bones has increased dramatically, the result of new cartilage generation between the bones. The outcome has been pain free walking for this 70-year-old female patient.




X-ray showing cartilage degeneration (left) and regeneration (right) using PEMF therapy.

Software package included with Curatron PC system


- All therapy parameters are under direct command and control of the program running on Windows/Mac.
- User-friendly database for demographic patient data including treatment follow up.
- Therapy settings can be selected from a database containing an extensive list of almost one hundred ready for use therapy protocols.
- In addition to the “ready for use treatment protocols,” all the therapy parameters can be set manually and saved in the Custom Disorders database and then run automatically.
- Free notes editing.
- Full PEMF report printing, including treatment history data, treatment time and accumulated treatment time.
- Automatic USB communication port detection for very easy installation without basic computer knowledge.

Easy to Use Software for PC or MAC platform

Record Patient Data and treatments. ▼



Choose a disorder from the list of pre-programmed settings. ►



Create a custom disorder with the software, and save your settings for future use. ►



Monitor Patient Progress ▼



Produce Patient Reports ▼

PEMF Software Patient Report

Patient Details

Name: James Bond

Birth date: 007

Notes:

Suffering from multiple gun wounds.

Diagnosed with Complex Regional Pain Syndrome resulting from severe trauma at lower left leg.

Created treatment protocol consisting of a series of frequencies with increasing intensity levels followed by high intensity local treatment with 8 Hz for 3 minutes.

Follow up after 3 weeks and repeated pain scale test (currently 8).

Date	Disorder	Duration (HH:MM:SS)
18-Sep-14 17:03:30	CRPS	0:00:52
18-Sep-14 17:00:15	CRPS	0:02:39

Total treatment time: 0:03:31

Applicators for Different Purposes:

The 2000PC works with a variety of applicators options including:

- Therapy mattress (26"x71")
- Therapy pad (20"x28")
- High energy coil (8"x8")
- Very high energy coil (8"x8")
- Hand-held probe



Very High Energy Coil for concentrated treatment of specific areas.



An optional foldable travel mat is available for use with the 2000PC and can be easily transported with the device.



Curatronic Ltd.

2000PC Specifications

2000PC System	Optional Applicators	Maximum Intensity
Energy Mattress	26"x71" (65cm x 180cm)	15,000 microTesla
Energy Pad	20"x28" (50cm x 70cm)	40,000 microTesla
High Energy Coil (A)	8"x8" (20cm x 20cm)	70,000 microTesla
Very High Energy Coil (B)	8"x8" (20cm x 20cm)	100,000 microTesla
Probe	6" x 1" (15cm x 2.5cm)	10,000 microTesla

Each organ in the body has cells which operate at different frequencies. Curatronic® conducts research into specific cell frequency and has preprogrammed its devices to stimulate cells at their optimum frequency related to specific disorders. The following is a list disorders and syndromes which the Curatronic® PC has preprogrammed settings for PEMF treatment on its devices:

Programmed Disorders:

- Acne
- Alzheimer's
- Ankylosing Spondylitis
- Anxiety
- Arthritis
- Atherosclerosis
- Asthma
- Athletic fatigue
- Backache
- Brachial neuralgia
- Bronchitis
- Burns
- Bursitis
- Carpal tunnel syndrome
- Cervical arthritis
- Cervical myalgic headaches
- Cervico brachialgia
- Coxarthrosis
- Coxitis
- Complex regional pain syndrome
- Contusions

- Decubitus
- Delayed fracture healing
- Dental treatment
- Dermatitis
- Dislocations
- Epicondylitis
- Facial neuralgia
- Fibromyalgia
- Fracture healing
- Gonarthrosis
- Herpes
- Hallux valgus bursitis
- Heart
- Intercostal neuralgia
- Intramuscular fibrosis
- Internal injuries
- Irritable bowel Syndrome
- Kidney resonance
- Liver resonance
- Low back pain
- Metabolic deficiency
- Migraine
- Migrainous neuralgia
- Morbus bechterew
- Morbus Sudeck
- Neuritis
- Neurosclerosis
- Neurosis
- Neurovegetative somatic symptoms
- Oedema
- Osteoarthritis
- Osteoporosis
- Osteomyelitis
- Pain therapy
- Pancreatitis
- Periarthritis
- Phlebitis
- Phobia
- Post radiotherapy treatment
- Prostatitis
- Prostate hypertrophy
- Psuedo arthrosis
- Psoriasis
- Relaxation
- Relaxation EEG waves
- Retinitis
- Rheumatoid arthritis
- Sacral arthritis

- Scars
- Sedation
- Sinusitis
- Sleep disturbances
- Spastic colon
- Stomach relaxation
- Stress
- Stroke recovery
- Tendinitis
- Thigh neuralgia
- Tonsillitis
- Torn muscles
- Trigeminus neuralgia
- Varicose veins
- Vitalization
- Wellness
- Wound healing

Brain Disorders:

- Alzheimer's
- Depression
- Neuropathic pain
- Parkinson
- Post traumatic stress
- Smoking cessation

Curatronic is the leading manufacturer of PEMF treatment systems which has proven itself as the standard to meet in both performance and long lasting effects.

www.electro-magnetic-therapy.com
info@curatronic.com



Curatronic Ltd.

www.curatronic.com www.curatronic.blogspot.com
www.pemft.com www.facebook.com/PEMFT
www.curavet.com www.youtube.com/user/Curatron/videos



Disclaimer: Pulsed electromagnetic field therapy (PEMF) has been approved by health authorities as medical therapy for human applications. In the European Community the Curatron devices are certified as medical devices according to the Medical Device Directive 93/42/EEC. The devices are manufactured according to ISO 13485 for Good Manufacturing Practice. In the United States of America the Curatron devices have not been cleared by the FDA. The information provided on this website is not medical advice. Health Canada has not sanctioned PEMF devices, however, PEMF technology is used by hospitals throughout Canada for certain applications. The products marketed on the website and in the brochures are not intended to diagnose, treat, cure or prevent any disease. Readers are encouraged to obtain the help, services and recommendations of doctors and other licensed medical practitioners.