



Curatronic Ltd.

# For Health Care Practitioners or Home Care Therapy

**Superior Power**



# Curatronic 2000XP<sup>SE</sup> Special Edition



# Curatron 2000XP<sup>SE</sup>

The Curatron 2000XP is the most popular Curatron model used by health care practitioners and home care therapy for informed individuals.

It is based on state-of-the-art technology generating Frequency Modulated Pulsed ElectroMagnetic Fields (PEMF). These fields are completely safe and cause no side effects. Curatron devices have undergone multiple scientific clinical trials worldwide. Tens of thousands of patients have been treated successfully with these systems. A high success percentage is reported, creating a dramatic improvement in the quality of life.

The Curatron 2000XP system increases peripheral oxygen pressure, activating and regenerating cells. Optimal oxygenation in the body is obtained by this professional, high-energy system, resulting in maximum therapeutic effect.



The Curatron 2000XP model excels in rapid pain reduction and has proven to be enormously beneficial for chronic pain therapy. People suffering from arthritis, fibromyalgia and back pain report almost instantaneous pain relief!

## Bone and Cartilage Regeneration

Improved calcium transport increases the absorption of calcium in bones and improves the quality of cartilage in joints, decreasing pain dramatically. Substantial increase in bone density is obtained in cases of osteoporosis, which in turn decreases fracture risk. Bone growth is stimulated causing fractured bones to heal properly.



*X-ray showing cartilage degeneration (left) and regeneration (right) using PEMF therapy.*

On the X-ray on the left, the distance between the bones is almost non-existent. After PEMF treatment, the X-ray on the right shows the distance between the bones has increased dramatically, the result of new cartilage generation between the bones. The outcome has been pain free walking for this 70-year-old female patient.

## Real Therapeutic Effect

The Curatron 2000XP is the only economically priced computerized therapy device which uses frequency modulated pulsed electromagnetic field technology. It generates sufficient energy in combination with rapid speed of induction to obtain a measurable therapeutic effect. 10 Different therapy programs are readily available, allowing for ideal individual automatic parameter settings for the listed disorders.

The pulsed electromagnetic fields are applied to the body by means of multiple flat coils located inside a therapy pad. The coil is placed at the area to be treated and is hooked up to the Curatron unit via a flexible cable.

## Easy to Use

The controls of the Curatron system are very user-friendly and simple to operate. Select a program number and press the start button and the complete therapy session will be performed automatically, under full control of the built-in computer. After the therapy time is over, the system switches to stand-by mode ensuring user safety while relaxing or sleeping.

## Special Edition

The 2000XP Special Edition is upgraded to deliver the power of the 2000PC. It boosts the output of the mattress and coils to maximize the healing potential of the unit.

# General Applications

### Wellness

- Anti stress
- Improve fitness & performance
- Anti aging
- Sleep improvement
- Faster recovery after injuries
- Metabolic stimulation

### Therapy

- Pain relief
- Arthritis
- Fibromyalgia
- Osteoporosis
- Fracture healing
- Sports injuries
- Rheumatic pain
- Tendinitis
- Low back pain
- Neuralgia
- Wound healing



*Very High Energy Coil for concentrated treatment of specific areas.*

## Applicators for Different Purposes:

The 2000XP<sup>SE</sup> works with a variety of applicators options including:

- Therapy mattress (26"x71")
- Therapy pad (20"x28")
- High energy coil (8"x8")
- Very high energy coil (8"x8")
- Hand held probe

## Take your Curatron on the Road



An optional foldable travel mat is available for use with the 2000XP and can be easily transported with the device.



Curatronic Ltd.

## 2000XP<sup>SE</sup> Specifications

2000XP <sup>SE</sup> System	Optional Applicators	Maximum Intensity
Energy Mattress	26"x71" (65cm x 180cm)	15,000 microTesla
Energy Pad	20"x28" (50cm x 70cm)	40,000 microTesla
High Energy Coil (A)	8"x8" (20cm x 20cm)	70,000 microTesla
Very High Energy Coil (B)	8"x8" (20cm x 20cm)	100,000 microTesla
Probe	6" x 1" (15cm x 2.5cm)	10,000 microTesla

# Enjoy life to the fullest with PEMF!



Health and Wellness are two of the most important things we can invest in. Without these every other aspect of our lives is diminished. People should be able to enjoy every day whether at work, school, or spending time with loved ones. However, each of us will face medical challenges in life, during which PEMF can help significantly.

Curatron's PEMF devices are designed to stimulate wellness by ensuring maximum blood flow, oxygenation, nutritional uptake, and detoxification. This allows the body to perform the tasks it was designed for. The goal of the manufacturer and the distributor of Curatron systems is to maximize health ensuring a superior quality of life to be experienced.

When disease and sickness occur,

PEMF energy medicine can play a major supporting role, assisting the body's built in defense and repair mechanisms to perform with optimal efficiency.

Curatron has developed its series of devices to deliver effective doses of energy medicine. It has worked with doctors to build program settings designed to address specific conditions. In some cases the affected organs are targetted, energizing the body to speed recovery and repair.

A partner to a healthy lifestyle including regular exercise, and a responsible diet, Curatron's PEMF systems are a prudent investment for those who want to get the most out of the gift of life.



Curatronic is the leading manufacturer of PEMF treatment systems which has proven itself as the standard to meet in both performance and long lasting effects.

[www.electro-magnetic-therapy.com](http://www.electro-magnetic-therapy.com)

[info@curatronic.com](mailto:info@curatronic.com)



## Curatronic Ltd.

[www.curatronic.com](http://www.curatronic.com) [www.curatronic.blogspot.com](http://www.curatronic.blogspot.com)

[www.pemft.com](http://www.pemft.com) [www.facebook.com/PEMFT](http://www.facebook.com/PEMFT)

[www.curavet.com](http://www.curavet.com) [www.youtube.com/user/Curatron/videos](http://www.youtube.com/user/Curatron/videos)



Disclaimer: Pulsed electromagnetic field therapy (PEMF) has been approved by health authorities as medical therapy for human applications. In the European Community the Curatron devices are certified as medical devices according to the Medical Device Directive 93/42/EEC. The devices are manufactured according to ISO 13485 for Good Manufacturing Practice. In the United States of America the Curatron devices have not been cleared by the FDA. The information provided on this website is not medical advice. Health Canada has not sanctioned PEMF devices, however, PEMF technology is used by hospitals throughout Canada for certain applications. The products marketed on the website and in the brochures are not intended to diagnose, treat, cure or prevent any disease. Readers are encouraged to obtain the help, services and recommendations of doctors and other licensed medical practitioners.